



CALIBRATION SCREEN TRAY VERIFICATION





Based upon the data below, Gilson Company, Inc. **VERIFIES** that this screen tray meets the requirements set forth by ASTM in Standard Specification for Woven Wire Test Sieve Cloth and Test Sieves Designation: E11-16, for **Inspection Test Sieves**.

The data for this screen tray was derived using a View Voyager 12X12 optical imaging device, following internal work instruction # **WI-750-006**. Units of measurement are recorded in millimeters. The View Voyager 12X12 machine, serial # V121203060361 was last **calibrated on 10/6/2016** with calibration **due on 10/1/2017**. This optical calibration was performed using calibration standard, serial number: 89947-001 and is traceable to N.I.S.T by Test Report No. 683/284711-14.

- New certification
 Re-certification
 Sieve Visually Inspected
 Sieve cleaned

Wire Diameter

Opening Size



Sharon McLaughlin

11-10-16

Helen Santelle

11-10-16

Laboratory Technician

Date

Reviewed By

Date



- P 740.548.7298 / 800.444.1508
- F 740.548.5314 / 800.255.5314
- W www.globalgilson.com



VERIFICATION INFORMATION

This Screen Tray Frame with Wire Cloth assembly has been verified as a single unit to comply with the requirements of the specification(s) noted. Matching serial numbers traceable to verification documents are attached to both the screen tray and wire cloth. To remain in compliance, the entire assembly (both the wire cloth and tray frame) must bear the identical serial number. Replacement of wire cloth or tray frame will void all verification per noted specifications. If replacement wire cloth is required, please contact Gilson Customer Service. Gilson will re-verify the completed assembly before returning.

Gilson Company, Inc.
P.O. Box 200
Lewis Center, Ohio 43035-0200
Phone: 800-444-1508 or 740-548-7298
Fax: 800-255-5314 or 740-548-5314
E-mail: customerservice@gilsonco.com
Website: www.globalgilson.com

Revised: 09/08/2011